

RFPL4B Antibody



PACO04359

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

IHC:1:100-1:300

Protein Background:

Gene ID:

RFPL4B

Uniprot

Q6ZWI9

Synonyms:

RFPL4B; RNF211; Ret finger protein-like 4B; RING finger protein 211

Immunogen:

Synthesized peptide derived from the Internal region of human RFPL4B.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

N/A

N/A