

Phospho-MSN/RDX/EZR (T558) Antibody



PACO04362

Product Information

Size:

50ug

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC, IF

Recommended dilutions:

WB:1:500-1:2000, IHC:1:100-1:300,
IF:1:200-1:1000

Protein Background:**Gene ID:**

MSN/RDX/EZR

Uniprot

P26038/P35241/P15311

Synonyms:

MSN; Moesin; Membrane-organizing extension spike protein; RDX; Radixin; EZR; VIL2; Ezrin; Cytovillin; Villin-2; p81

Immunogen:

Synthesized peptide derived from human Moesin/Ezrin/Radixin around the phosphorylation site of T558.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

N/A

N/A