

Phospho-MYB (S532) Antibody



PACO04372

Product Information

Size:

50ug

Reactivity:

Human, Mouse

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

IHC:1:100-1:300

Protein Background:**Gene ID:**

MYB

Uniprot

P10242

Synonyms:

MYB; Transcriptional activator Myb; Proto-oncogene c-Myb

Immunogen:

Synthesized peptide derived from human c-Myb around the phosphorylation site of S532.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

N/A

N/A