

# Phospho-MYC (T358) Antibody



PACO04376

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human, Mouse, Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, IHC, IP

**Recommended dilutions:**

IHC:1:100-1:300, IP:2-5ug

**Protein Background:****Gene ID:**

MYC

**Uniprot**

P01106

**Synonyms:**

MYC; BHLHE39; Myc proto-oncogene protein; Class E basic helix-loop-helix protein 39; bHLHe39; Proto-oncogene c-Myc; Transcription factor p64

**Immunogen:**

Synthesized peptide derived from human c-Myc around the phosphorylation site of T358.

**Storage:**

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

## Product Images

---

N/A

N/A