

Phospho-PPP1R12A (T696) Antibody



PACO04391

Product Information

Size:

50ug

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC, IF

Recommended dilutions:

IHC:1:100-1:300, IF:1:200-1:1000

Protein Background:**Gene ID:**

PPP1R12A

Uniprot

O14974

Synonyms:

PPP1R12A; MBS; MYPT1; Protein phosphatase 1 regulatory subunit 12A; Myosin phosphatase-targeting subunit 1; Myosin phosphatase target subunit 1; Protein phosphatase myosin-binding subunit

Immunogen:

Synthesized peptide derived from human MYPT1 around the phosphorylation site of T696.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

N/A

N/A