

NKTR Antibody



PACO04446

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC, IF

Recommended dilutions:

IHC:1:100-1:300, IF:1:200-1:1000

Protein Background:

Gene ID:

NKTR

Uniprot

P30414

Synonyms:

NKTR; NK-tumor recognition protein; NK-TR protein; Natural-killer cells cyclophilin-related protein

Immunogen:

Synthesized peptide derived from the Internal region of human NK-TR.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

N/A

N/A