

NPY Antibody



PACO04456

Product Information

Size:

50ug

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

IHC:1:100-1:300, IHC:1:200-1:1000

Protein Background:**Gene ID:**

NPY

Uniprot

P01303

Synonyms:

NPY; Pro-neuropeptide Y

Immunogen:

Synthesized peptide derived from the C-terminal region of human NPY.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

N/A

N/A