

NPHS1 Antibody



PACO04459

Product Information

Size:

50ug

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

IHC:1:100-1:300

Protein Background:**Gene ID:**

NPHS1

Uniprot

O60500

Synonyms:

NPHS1; NPHN; Nephlin; Renal glomerulus-specific cell adhesion receptor

Immunogen:

Synthesized peptide derived from the C-terminal region of human Nephlin.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

N/A

N/A