

# ZBTB21 Antibody



PACO04481

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human, Mouse

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, IHC

**Recommended dilutions:**

IHC:1:100-1:300

**Protein Background:**

**Gene ID:**

ZBTB21

**Uniprot**

Q9ULJ3

**Synonyms:**

ZNF295; KIAA1227; ZBTB21; Zinc finger protein 295; Zinc finger and BTB domain-containing protein 21

**Immunogen:**

Synthesized peptide derived from the Internal region of human ZNF295.

**Storage:**

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

## Product Images

---

N/A

N/A