

OVOL1 Antibody



PACO04488

Product Information

Size:

50ug

Reactivity:

Human, Mouse

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

WB:1:500-1:2000, IHC:1:100-1:300

Protein Background:**Gene ID:**

OVOL1

Uniprot

O14753

Synonyms:

OVOL1; Putative transcription factor Ovo-like 1; hOvo1

Immunogen:

Synthesized peptide derived from the N-terminal region of human Ovol1.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

N/A

N/A