

TAS2R1 Antibody



PACO04507

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IF

Recommended dilutions:

WB:1:500-1:2000, IF:1:200-1:1000

Protein Background:

Gene ID:

TAS2R1

Uniprot

Q9NYW7

Synonyms:

TAS2R1; Taste receptor type 2 member 1; T2R1; Taste receptor family B member 7; TRB7

Immunogen:

Synthesized peptide derived from the N-terminal region of human T2R1.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

N/A

N/A