

# MRPS16 Antibody



PACO04516

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human, Mouse

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB, IHC

**Recommended dilutions:**

WB:1:500-1:2000, IHC:1:100-1:300

**Protein Background:**

**Gene ID:**

MRPS16

**Uniprot**

Q9Y3D3

**Synonyms:**

MRPS16; RPMS16; CGI-132; 28S ribosomal protein S16; mitochondrial; MRP-S16; S16mt

**Immunogen:**

Synthesized peptide derived from the Internal region of human MRP-S16.

**Storage:**

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

## Product Images

---

N/A

N/A