

CEND1 Antibody



PACO04528

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

IHC:1:100-1:300

Protein Background:**Gene ID:**

CEND1

Uniprot

Q8N111

Synonyms:

CEND1; BM88; Cell cycle exit and neuronal differentiation protein 1; BM88 antigen

Immunogen:

Synthesized peptide derived from the N-terminal region of human BM88.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

N/A

N/A