

# TMBIM4 Antibody



PACO04561

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human, Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB, IHC

**Recommended dilutions:**

WB:1:500-1:2000, IHC:1:100-1:300

**Protein Background:****Gene ID:**

TMBIM4

**Uniprot**

Q9HC24

**Synonyms:**

TMBIM4; GAAP; LFG4; CGI-119; Protein lifeguard 4; Golgi anti-apoptotic protein; Protein S1R; Transmembrane BAX inhibitor motif-containing protein 4; Z-protein

**Immunogen:**

Synthesized peptide derived from the Internal region of human GAAP.

**Storage:**

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

## Product Images

---

N/A

N/A