

Phospho-PDPK1 (S241) Antibody



PACO04565

Product Information

Size:

50ug

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

WB:1:500-1:2000, IHC:1:100-1:300

Protein Background:

Gene ID:

PDPK1

Uniprot

O15530

Synonyms:

PDPK1; PDK1; 3-phosphoinositide-dependent protein kinase 1; hPDK1

Immunogen:

Synthesized peptide derived from human PDK1 around the phosphorylation site of S241.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

N/A

N/A