

PHAX Antibody



PACO04576

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

WB:1:500-1:2000, IHC:1:200-1:1000

Protein Background:

Gene ID:

PHAX

Uniprot

Q9H814

Synonyms:

PHAX; RNUXA; Phosphorylated adapter RNA export protein; RNA U small nuclear RNA export adapter protein

Immunogen:

Synthesized peptide derived from the Internal region of human PHAX.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

N/A

N/A