

# PIGH Antibody



PACO04587

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human, Mouse

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, IHC

**Recommended dilutions:**

IHC:1:100-1:300

**Protein Background:**

**Gene ID:**

PIGH

**Uniprot**

Q14442

**Synonyms:**

PIGH; Phosphatidylinositol N-acetylglucosaminyltransferase subunit H;  
Phosphatidylinositol-glycan biosynthesis class H protein; PIG-H

**Immunogen:**

Synthesized peptide derived from the C-terminal region of human PIG-H.

**Storage:**

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

## Product Images

---

N/A

N/A