

PIM1 Antibody



PACO04589

Product Information

Size:

50ug

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

WB:1:500-1:2000

Protein Background:**Gene ID:**

PIM1

Uniprot

P11309

Synonyms:

PIM1; Serine/threonine-protein kinase pim-1

Immunogen:

Synthesized peptide derived from human Pim-1 around the non-phosphorylation site of Y309.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

N/A

N/A