

Phospho-PLN (S16/T17) Antibody



PACO04615

Product Information

Size:

50ug

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IF

Recommended dilutions:

IF:1:200-1:1000

Protein Background:

Gene ID:

PLN

Uniprot

P26678

Synonyms:

PLN; PLB; Cardiac phospholamban; PLB

Immunogen:

Synthesized peptide derived from human PLB around the phosphorylation site of S16/T17.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

N/A

N/A