

**Phospho-CSNK1G1/CSNK1G2/CSNK1G3 (Y263)
Antibody
PACO04623**



Product Information

Size:

50ug

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

IHC:1:100-1:300

Protein Background:

Gene ID:

CSNK1G1/CSNK1G2/CSNK1G3

Uniprot

Q9HCP0/P78368/Q9Y6M4

Synonyms:

CSNK1G1; Casein kinase I isoform gamma-1; CKI-gamma 1; CSNK1G2; CK1G2; Casein kinase I isoform gamma-2; CKI-gamma 2; CSNK1G3; Casein kinase I isoform gamma-3; CKI-gamma 3

Immunogen:

Synthesized peptide derived from human Casein Kinase I gamma1/2/3 around the phosphorylation site of Y263.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

N/A

N/A