

CLIC6 Antibody



PACO04627

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

WB:1:500-1:2000

Protein Background:

Gene ID:

CLIC6

Uniprot

Q96NY7

Synonyms:

CLIC6; CLIC1L; Chloride intracellular channel protein 6; Parchorin

Immunogen:

Synthesized peptide derived from the Internal region of human CLIC6.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

N/A

N/A