

CYTL1 Antibody



PACO04636

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

WB:1:500-1:2000, IHC:1:100-1:300

Protein Background:**Gene ID:**

CYTL1

Uniprot

Q9NRR1

Synonyms:

CYTL1; C4orf4; Cytokine-like protein 1; Protein C17

Immunogen:

Synthesized peptide derived from the Internal region of human CYTL1.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

N/A

N/A