

PRPF39 Antibody



PACO04670

Product Information

Size:

50ug

Reactivity:

Human, Mouse

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

WB:1:500-1:2000, IHC:1:100-1:300

Protein Background:**Gene ID:**

PRPF39

Uniprot

Q86UA1

Synonyms:

PRPF39; Pre-mRNA-processing factor 39; PRP39 homolog

Immunogen:

Synthesized peptide derived from the Internal region of human PRPF39.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

N/A

N/A