

# Phospho-PPP1R2 (S44) Antibody



PACO04673

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human, Mouse, Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB, IHC

**Recommended dilutions:**

WB:1:500-1:2000, IHC:1:100-1:300

**Protein Background:****Gene ID:**

PPP1R2

**Uniprot**

P41236

**Synonyms:**

PPP1R2; IPP2; Protein phosphatase inhibitor 2; IPP-2

**Immunogen:**

Synthesized peptide derived from human IPP-2 around the phosphorylation site of S44.

**Storage:**

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

## Product Images

---

N/A

N/A