

PPP1R2 Antibody



PACO04676

Product Information

Size:

50ug

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

WB:1:500-1:2000, IHC:1:100-1:300

Protein Background:**Gene ID:**

PPP1R2

Uniprot

P41236

Synonyms:

PPP1R2; IPP2; Protein phosphatase inhibitor 2; IPP-2

Immunogen:

Synthesized peptide derived from human IPP-2 around the non-phosphorylation site of S120/S121.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

N/A

N/A