

# PLCXD1 Antibody



PACO04695

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, IHC

**Recommended dilutions:**

IHC:1:100-1:300

**Protein Background:**

**Gene ID:**

PLCXD1

**Uniprot**

Q9NUJ7

**Synonyms:**

PLCXD1; PI-PLC X domain-containing protein 1

**Immunogen:**

Synthesized peptide derived from the Internal region of human PLC-XD1.

**Storage:**

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

## Product Images

---

N/A

N/A