

PPT1 Antibody



PACO04697

Product Information

Size:

50ug

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

IHC:1:100-1:300

Protein Background:**Gene ID:**

PPT1

Uniprot

P50897

Synonyms:

PPT1; PPT; Palmitoyl-protein thioesterase 1; PPT-1; Palmitoyl-protein hydrolase 1

Immunogen:

Synthesized peptide derived from the N-terminal region of human CLN1.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

N/A

N/A