

CSGALNACT2 Antibody



PACO04700

Product Information

Size:

50ug

Reactivity:

Human, Mouse

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

IHC:1:100-1:300

Protein Background:

Gene ID:

CSGALNACT2

Uniprot

Q8N6G5

Synonyms:

CSGALNACT2; CHGN2; GALNACT2; Chondroitin sulfate N-acetylgalactosaminyltransferase 2; Chondroitin beta-1; 4-N-acetylgalactosaminyltransferase 2; Beta4GalNAcT-2; GalNAcT-2

Immunogen:

Synthesized peptide derived from the Internal region of human CsGalNAcT-2.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

N/A

N/A