

ZNF446 Antibody



PACO04708

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

WB:1:500-1:2000

Protein Background:**Gene ID:**

ZNF446

Uniprot

Q9NWS9

Synonyms:

ZNF446; ZKSCAN20; Zinc finger protein 446; Zinc finger protein with KRAB and SCAN domains 20

Immunogen:

Synthesized peptide derived from the Internal region of human ZNF446.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

N/A

N/A