

PRKACB Antibody



PACO04710

Product Information

Size:

50ug

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

WB:1:500-1:2000, IHC:1:100-1:300

Protein Background:**Gene ID:**

PRKACB

Uniprot

P22694

Synonyms:

PRKACB; cAMP-dependent protein kinase catalytic subunit beta; PKA C-beta

Immunogen:

Synthesized peptide derived from the C-terminal region of human PKABeta cat.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

N/A

N/A