

PACO04711

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human, Mouse, Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, IHC, IF

**Recommended dilutions:**

IHC:1:100-1:300, IF:1:200-1:1000

**Protein Background:****Gene ID:**

PKIA

**Uniprot**

P61925

**Synonyms:**

PKIA; PRKACN1; cAMP-dependent protein kinase inhibitor alpha; PKI-alpha; cAMP-dependent protein kinase inhibitor; muscle/brain isoform

**Immunogen:**

Synthesized peptide derived from the N-terminal region of human Pki alpha.

**Storage:**

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

## Product Images

---

N/A

N/A