

Phospho-PRKAR2B (S113) Antibody



PACO04721

Product Information

Size:

50ug

Reactivity:

Human, Mouse, Rat, Monkey

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

WB:1:500-1:2000, IHC:1:100-1:300

Protein Background:**Gene ID:**

PRKAR2B

Uniprot

P31323

Synonyms:

PRKAR2B; cAMP-dependent protein kinase type II-beta regulatory subunit

Immunogen:

Synthesized peptide derived from human PKA IIBeta reg around the phosphorylation site of S113.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

N/A

N/A