

# ZNF839 Antibody



PACO04722

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, IHC

**Recommended dilutions:**

IHC:1:100-1:300

**Protein Background:**

**Gene ID:**

ZNF839

**Uniprot**

A8K0R7

**Synonyms:**

ZNF839; C14orf131; Zinc finger protein 839; Renal carcinoma antigen NY-REN-50

**Immunogen:**

Synthesized peptide derived from the Internal region of human NY-REN-50.

**Storage:**

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

## Product Images

---

N/A

N/A