

CSGALNACT1 Antibody



PACO04727

Product Information

Size:

50ug

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

WB:1:500-1:2000, IHC:1:100-1:300

Protein Background:**Gene ID:**

CSGALNACT1

Uniprot

Q8TDX6

Synonyms:

CSGALNACT1; CHGN; GALNACT1; Chondroitin sulfate N-acetylgalactosaminyltransferase 1; CsGalNACT-1; Chondroitin beta-1; 4-N-acetylgalactosaminyltransferase 1; Beta4GalNACT-1

Immunogen:

Synthesized peptide derived from the Internal region of human Beta-1,4-GalNAc-T.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

N/A

N/A