

# Phospho-PRKCQ (S676) Antibody



PACO04739

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human, Mouse, Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB, IHC, IF

**Recommended dilutions:**

WB:1:500-1:2000, IHC:1:100-1:300,  
IF:1:200-1:1000

**Protein Background:**

**Gene ID:**

PRKCQ

**Uniprot**

Q04759

**Synonyms:**

PRKCQ; PRKCT; Protein kinase C theta type; nPKC-theta

**Immunogen:**

Synthesized peptide derived from human PKC Theta around the phosphorylation site of S676.

**Storage:**

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

## Product Images

---

N/A

N/A