

Phospho-PRKCQ (T538) Antibody



PACO04740

Product Information

Size:

50ug

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

WB:1:500-1:2000, IHC:1:100-1:300

Protein Background:**Gene ID:**

PRKCQ

Uniprot

Q04759

Synonyms:

PRKCQ; PRKCT; Protein kinase C theta type; nPKC-theta

Immunogen:

Synthesized peptide derived from human PKC Theta around the phosphorylation site of T538.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

N/A

N/A