

Phospho-PRKDC (T2638) Antibody



PACO04748

Product Information

Size:

50ug

Reactivity:

Human, Mouse

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IF

Recommended dilutions:

IF:1:200-1:1000

Protein Background:**Gene ID:**

PRKDC

Uniprot

P78527

Synonyms:

PRKDC; HYRC; HYRC1; DNA-dependent protein kinase catalytic subunit; DNA-PK catalytic subunit; DNA-PKcs; DNPk1; p460

Immunogen:

Synthesized peptide derived from human DNA-PKCS around the phosphorylation site of T2638.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

N/A

N/A