

NXF3 Antibody



PACO04760

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

WB:1:500-1:2000, IHC:1:100-1:300

Protein Background:

Gene ID:

NXF3

Uniprot

Q9H4D5

Synonyms:

NXF3; TAPL3; Nuclear RNA export factor 3; TAP-like protein 3; TAPL-3

Immunogen:

Synthesized peptide derived from the Internal region of human NXF3.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

N/A

N/A