

FMN2 Antibody



PACO04785

Product Information

Size:

50ug

Reactivity:

Human, Mouse

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

IHC:1:100-1:300

Protein Background:

Gene ID:

FMN2

Uniprot

Q9NZ56

Synonyms:

FMN2; Formin-2

Immunogen:

Synthesized peptide derived from the C-terminal region of human Formin 2.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

N/A

N/A