

PTGIS Antibody



PACO04809

Product Information

Size:

50ug

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

WB:1:500-1:2000

Protein Background:

Gene ID:

PTGIS

Uniprot

Q16647

Synonyms:

PTGIS; CYP8; CYP8A1; Prostacyclin synthase; Prostaglandin I2 synthase

Immunogen:

Synthesized peptide derived from the C-terminal region of human PGI2 synthase.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

N/A

N/A