

BAI3 Antibody



PACO04834

Product Information

Size:

50ug

Reactivity:

Human, Mouse

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC, IF

Recommended dilutions:

IHC:1:100-1:300, IF:1:200-1:1000

Protein Background:**Gene ID:**

ADGRB3

Uniprot

O60242

Synonyms:

BAI3; KIAA0550; Brain-specific angiogenesis inhibitor 3

Immunogen:

Synthesized peptide derived from the N-terminal region of human BAI-3.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

N/A

N/A