

# Phospho-RAD52 (Y104) Antibody



PACO04871

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human, Mouse, Monkey

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB, IHC

**Recommended dilutions:**

WB:1:500-1:2000, IHC:1:100-1:300

**Protein Background:**

**Gene ID:**

RAD52

**Uniprot**

P43351

**Synonyms:**

RAD52; DNA repair protein RAD52 homolog

**Immunogen:**

Synthesized peptide derived from human Rad52 around the phosphorylation site of Y104.

**Storage:**

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

## Product Images

---

N/A

N/A