

RB1 Antibody



PACO04892

Product Information

Size:

50ug

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC, IF

Recommended dilutions:

IHC:1:100-1:300, IF:1:200-1:1000

Protein Background:**Gene ID:**

RB1

Uniprot

P06400

Synonyms:

RB1; Retinoblastoma-associated protein; p105-Rb; pRb; Rb; pp110

Immunogen:

Synthesized peptide derived from human Rb around the non-phosphorylation site of T826.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

N/A

N/A