

Phospho-BCL2 (T74) Antibody



PACO04903

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

IHC:1:100-1:300

Protein Background:

Gene ID:

BCL2

Uniprot

P10415

Synonyms:

BCL2; Apoptosis regulator Bcl-2

Immunogen:

Synthesized peptide derived from human Bcl-2 around the phosphorylation site of T74.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

N/A

N/A