

Phospho-RHO (S334) Antibody



PACO04923

Product Information

Size:

50ug

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

IHC:1:100-1:300

Protein Background:

Gene ID:

RHO

Uniprot

P08100

Synonyms:

RHO; OPN2; Rhodopsin; Opsin-2

Immunogen:

Synthesized peptide derived from human Rhodopsin around the phosphorylation site of S334.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

N/A

N/A