

RPL11 Antibody



PACO04941

Product Information

Size:

50ug

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

WB:1:500-1:2000

Protein Background:**Gene ID:**

RPL11

Uniprot

P62913

Synonyms:

RPL11; 60S ribosomal protein L11; CLL-associated antigen KW-12

Immunogen:

Synthesized peptide derived from the C-terminal region of human Ribosomal Protein L11.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

N/A

N/A