

Phospho-RPS6 (S240) Antibody



PACO04951

Product Information

Size:

50ug

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC, IF

Recommended dilutions:

WB:1:500-1:2000, IHC:1:100-1:300,
IF:1:200-1:1000

Protein Background:**Gene ID:**

RPS6

Uniprot

P62753

Synonyms:

RPS6; OK/SW-cl.2; 40S ribosomal protein S6; Phosphoprotein NP33

Immunogen:

Synthesized peptide derived from human Ribosomal Protein S6 around the phosphorylation site of S240.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

N/A

N/A