

# Phospho-RPS6KA1/RPS6KA3/RPS6KA2/RPS6KA6 (S221/227/S218/232) Antibody

PACO04956



---

## Product Information

**Size:**

50ug

**Reactivity:**

Human, Mouse

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB

**Recommended dilutions:**

WB:1:500-1:2000

**Protein Background:****Gene ID:**

RPS6KA1/RPS6KA3/RPS6KA2/RPS6KA6

**Uniprot**

Q15418/P51812/Q15349/Q9UK32

**Synonyms:**

RPS6KA1; MAPKAPK1A; RSK1; Ribosomal protein S6 kinase alpha-1; S6K-alpha-1; 90 kDa ribosomal protein S6 kinase 1; p90-RSK 1; p90RSK1; p90S6K; MAP kinase-activated protein kinase 1a; MAPK-activated protein kinase 1a; MAPKAP kinase 1a; MAPKAP

**Immunogen:**

Synthesized peptide derived from human Rsk-1/2/3/4 around the phosphorylation site of S221/227/S218/232.

**Storage:**

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

## Product Images

---

N/A

N/A