

RPS6KA1/RPS6KA3/RPS6KA2/RPS6KA6 Antibody



PACO04957

Product Information

Size:

50ug

Reactivity:

Human, Mouse

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

WB:1:500-1:2000

Protein Background:**Gene ID:**

RPS6KA1/RPS6KA3/RPS6KA2/RPS6KA6

Uniprot

Q15418/P51812/Q15349/Q9UK32

Synonyms:

RPS6KA1; MAPKAPK1A; RSK1; Ribosomal protein S6 kinase alpha-1; S6K-alpha-1; 90 kDa ribosomal protein S6 kinase 1; p90-RSK 1; p90RSK1; p90S6K; MAP kinase-activated protein kinase 1a; MAPK-activated protein kinase 1a; MAPKAP kinase 1a; MAPKAP

Immunogen:

Synthesized peptide derived from human Rsk-1/2/3/4 around the non-phosphorylation site of S221/227/S218/232.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

N/A

N/A