

# RPS18 Antibody



PACO04968

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human, Mouse, Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB

**Recommended dilutions:**

WB:1:500-1:2000

**Protein Background:****Gene ID:**

RPS18

**Uniprot**

P62269

**Synonyms:**

RPS18; D6S218E; 40S ribosomal protein S18; Ke-3; Ke3

**Immunogen:**

Synthesized peptide derived from the Internal region of human Ribosomal Protein S18.

**Storage:**

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

## Product Images

---

N/A

N/A